

V-set domain containing T-cell activation inhibitor 1

Catalog No: #AP89337



Package Size: #AP89337-1 50ug #AP89337-2 100ug #AP89337-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	V-set domain containing T-cell activation inhibitor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B7 homolog 4,B7-H4,Immune costimulatory protein B7-H4,Protein B7S1,T cell costimulatory molecule B7x,B7h4
Accession No.	Uniprot ID:Q7TSP5
Uniprot	Q7TSP5
GeneID	242122;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vtcn1

Note: This product is for in vitro research use only