

# WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1

Catalog No: #AP89338

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP89338-1 50ug #AP89338-2 100ug #AP89338-3 1mg

## Description

Product Name	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Growth and differentiation factor-associated serum protein 2,GASP-2,Gasp2
Accession No.	Uniprot ID:Q8R0S6
Uniprot	Q8R0S6
GeneID	215001;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Wfikkn1

Note: This product is for in vitro research use only