

Disintegrin and metalloproteinase domain-containing protein 29



Catalog No: #AP89339

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89339-1 50ug #AP89339-2 100ug #AP89339-3 1mg

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 29
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM 29,Cancer/testis antigen 73,CT73
Accession No.	Uniprot ID:Q9UKF5
Uniprot	Q9UKF5
GeneID	11086;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ADAM29

Note: This product is for in vitro research use only