

Xanthine dehydrogenase/oxidase

Catalog No: #AP89340



Package Size: #AP89340-1 50ug #AP89340-2 100ug #AP89340-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Xanthine dehydrogenase/oxidase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Accession No. | Uniprot ID:P80457 |
| Uniprot | P80457 |
| GeneID | 280960; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:XDH

Note: This product is for in vitro research use only