

Adenylate cyclase type 10

Catalog No: #AP89341



Package Size: #AP89341-1 50ug #AP89341-2 100ug #AP89341-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenylate cyclase type 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Germ cell soluble adenylyl cyclase,sAC,Testicular soluble adenylyl cyclase,Sac,Sacy
Accession No.	Uniprot ID:Q8C0T9
Uniprot	Q8C0T9
GeneID	271639;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adcy10

Note: This product is for in vitro research use only