

## ADP/ATP translocase 3

Catalog No: #AP89342



Package Size: #AP89342-1 50ug #AP89342-2 100ug #AP89342-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ADP/ATP translocase 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADP,ATP carrier protein 3,ADP,ATP carrier protein, isoform T2,ANT 2,Adenine nucleotide translocator 3,ANT 3,Solute carrier family 25 member 6,ANT3
Accession No.	Uniprot ID:P32007
Uniprot	P32007
GeneID	282480;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC25A6

Note: This product is for in vitro research use only