

Angiotensin-converting enzyme 2

Catalog No: #AP89345

Package Size: #AP89345-1 50ug #AP89345-2 100ug #AP89345-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Angiotensin-converting enzyme 2, ACE2
Accession No.	Uniprot ID:Q2TNK5
Uniprot	Q2TNK5
GeneID	397628;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ANGPTL4

Note: This product is for in vitro research use only