

Aprataxin and PNK-like factor

Catalog No: #AP89348



Package Size: #AP89348-1 50ug #AP89348-2 100ug #AP89348-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aprataxin and PNK-like factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Apurinic-apyrimidinic endonuclease APLF
Accession No.	Uniprot ID:A0JNH9
Uniprot	A0JNH9
GeneID	507989;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APLF

Note: This product is for in vitro research use only