

Apolipoprotein C-IV

Catalog No: #AP89350



Package Size: #AP89350-1 50ug #AP89350-2 100ug #AP89350-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein C-IV
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Apo-CIV, Apolipoprotein C4
Accession No.	Uniprot ID:Q3SYR5
Uniprot	Q3SYR5
GeneID	618041;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APOC4

Note: This product is for in vitro research use only