

## S-arrestin

Catalog No: #AP89352



Package Size: #AP89352-1 50ug #AP89352-2 100ug #AP89352-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	S-arrestin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	48 kDa protein, Retinal S-antigen, S-AG, Rod photoreceptor arrestin
Accession No.	Uniprot ID:P20443
Uniprot	P20443
GeneID	20215;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Sag

Note: This product is for in vitro research use only