

Cytoplasmic tyrosine-protein kinase BMX

Catalog No: #AP89356



Package Size: #AP89356-1 50ug #AP89356-2 100ug #AP89356-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytoplasmic tyrosine-protein kinase BMX
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Bone marrow tyrosine kinase gene in chromosome X protein,Epithelial and endothelial tyrosine kinase,ETK,NTK38
Accession No.	Uniprot ID:P51813
Uniprot	P51813
GeneID	660;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BMX

Note: This product is for in vitro research use only