

C-C motif chemokine 4

Catalog No: #AP89363

Package Size: #AP89363-1 50ug #AP89363-2 100ug #AP89363-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Macrophage inflammatory protein 1-beta,MIP-1-beta,Small-inducible cytokine A4,MIP1B,SCY4A
Accession No.	Uniprot ID:Q711P4
Uniprot	Q711P4
GeneID	396668;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL4

Note: This product is for in vitro research use only