

Casein kinase II subunit alpha

Catalog No: #AP89375



Package Size: #AP89375-1 50ug #AP89375-2 100ug #AP89375-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Casein kinase II subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CK II alpha,CK2A1
Accession No.	Uniprot ID:P68399
Uniprot	P68399
GeneID	282419;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CSNK2A1

Note: This product is for in vitro research use only