

Glutamate decarboxylase 1

Catalog No: #AP89378



Package Size: #AP89378-1 50ug #AP89378-2 100ug #AP89378-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate decarboxylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	67 kDa glutamic acid decarboxylase, GAD-67, Glutamate decarboxylase 67 kDa isoform, Gad67
Accession No.	Uniprot ID:P48318
Uniprot	P48318
GeneID	14415;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gad1

Note: This product is for in vitro research use only