

DNA damage-inducible transcript 3 protein

Catalog No: #AP89380



Package Size: #AP89380-1 50ug #AP89380-2 100ug #AP89380-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA damage-inducible transcript 3 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	DDIT-3,C/EBP zeta,C/EBP-homologous protein,CHOP,C/EBP-homologous protein 10,CHOP-10,CCAAT/enhancer-binding protein homologous protein,Growth arrest and DNA-damage-inducible protein GADD153,CHOP,CHOP10,GADD153
Accession No.	Uniprot ID:Q0IIB6
Uniprot	Q0IIB6
GeneID	777788;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDIT3

Note: This product is for in vitro research use only