

Endoplasmic reticulum metallopeptidase 1

Catalog No: #AP89386



Package Size: #AP89386-1 50ug #AP89386-2 100ug #AP89386-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endoplasmic reticulum metallopeptidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Felix-ina,D19Wsu12e,Fxna,Kiaa1815
Accession No.	Uniprot ID:Q3UVK0
Uniprot	Q3UVK0
GeneID	226090;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ermp1

Note: This product is for in vitro research use only