

Fatty acid desaturase 1

Catalog No: #AP89389



Package Size: #AP89389-1 50ug #AP89389-2 100ug #AP89389-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid desaturase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Delta(5) fatty acid desaturase,D5D,FADSD5
Accession No.	Uniprot ID:O60427
Uniprot	O60427
GeneID	3992;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FADS1

Note: This product is for in vitro research use only