

Fibroblast growth factor 5

Catalog No: #AP89391



Package Size: #AP89391-1 50ug #AP89391-2 100ug #AP89391-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FGF-5,Heparin-binding growth factor 5,HBGF-5,Fgf-5
Accession No.	Uniprot ID:P48807
Uniprot	P48807
GeneID	60662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgf5

Note: This product is for in vitro research use only