

Peptidyl-prolyl cis-trans isomerase FKBP8

Catalog No: #AP89392



Package Size: #AP89392-1 50ug #AP89392-2 100ug #AP89392-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase FKBP8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PPlase FKBP8,FK506-binding protein 8,FKBP-8,Rotamase
Accession No.	Uniprot ID:Q3B7U9
Uniprot	Q3B7U9
GeneID	290652;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fkbp8

Note: This product is for in vitro research use only