

Receptor-type tyrosine-protein kinase FLT3

Catalog No: #AP89393



Package Size: #AP89393-1 50ug #AP89393-2 100ug #AP89393-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-type tyrosine-protein kinase FLT3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FL cytokine receptor,Fetal liver kinase-2,FLK-2,Fms-like tyrosine kinase 3,FLT-3,Stem cell tyrosine kinase 1,STK-1,CD135,CD135,FLK2,STK1
Accession No.	Uniprot ID:P36888
Uniprot	P36888
GeneID	2322;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FLT3

Note: This product is for in vitro research use only