

Forkhead box protein C2

Catalog No: #AP89394



Package Size: #AP89394-1 50ug #AP89394-2 100ug #AP89394-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Forkhead box protein C2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain factor 3,BF-3,Forkhead-related protein FKHL14,Mesenchyme fork head protein 1,MFH-1 protein,Transcription factor FKH-14,Fkh14,Fkh14,Mfh1
Accession No.	Uniprot ID:Q61850
Uniprot	Q61850
GeneID	14234;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Foxc2

Note: This product is for in vitro research use only