

Platelet glycoprotein V

Catalog No: #AP89399



Package Size: #AP89399-1 50ug #AP89399-2 100ug #AP89399-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Platelet glycoprotein V
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GPV,Glycoprotein 5,CD42d
Accession No.	Uniprot ID:P40197
Uniprot	P40197
GeneID	2814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GP5

Note: This product is for in vitro research use only