

Glutathione S-transferase P 1

Catalog No: #AP89401



Package Size: #AP89401-1 50ug #AP89401-2 100ug #AP89401-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutathione S-transferase P 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gst P1,GST YF-YF,GST class-pi,GST-piB,Preadipocyte growth factor,Gstpb
Accession No.	Uniprot ID:P19157
Uniprot	P19157
GeneID	14870;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gstp1

Note: This product is for in vitro research use only