

Porphobilinogen deaminase

Catalog No: #AP89403



Package Size: #AP89403-1 50ug #AP89403-2 100ug #AP89403-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Porphobilinogen deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PBG-D,Hydroxymethylbilane synthase,HMBS,Pre-uroporphyrinogen synthase,Uros1
Accession No.	Uniprot ID:P22907
Uniprot	P22907
GeneID	15288;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmbs

Note: This product is for in vitro research use only