

Beta-hexosaminidase subunit beta

Catalog No: #AP89404



Package Size: #AP89404-1 50ug #AP89404-2 100ug #AP89404-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-hexosaminidase subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-N-acetylhexosaminidase subunit beta, Hexosaminidase subunit B, N-acetyl-beta-glucosaminidase subunit beta
Accession No.	Uniprot ID:Q6AXR4
Uniprot	Q6AXR4
GeneID	294673;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hexb

Note: This product is for in vitro research use only