

Heparanase

Catalog No: #AP89405



Package Size: #AP89405-1 50ug #AP89405-2 100ug #AP89405-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heparanase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Endo-gluconidase,Hep
Accession No.	Uniprot ID:Q71RP1
Uniprot	Q71RP1
GeneID	64537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hpse

Note: This product is for in vitro research use only