

Histamine H4 receptor

Catalog No: #AP89406



Package Size: #AP89406-1 50ug #AP89406-2 100ug #AP89406-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histamine H4 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	H4R,AXOR35,G-protein coupled receptor 105,GPRv53,Pfi-013,SP9144,GPCR105
Accession No.	Uniprot ID:Q9H3N8
Uniprot	Q9H3N8
GeneID	59340;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HRH4

Note: This product is for in vitro research use only