

Heat shock-related 70 kDa protein 2

Catalog No: #AP89407



Package Size: #AP89407-1 50ug #AP89407-2 100ug #AP89407-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock-related 70 kDa protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Heat shock 70 kDa protein 2
Accession No.	Uniprot ID:P54652
Uniprot	P54652
GeneID	3306;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HSPA2

Note: This product is for in vitro research use only