

Interferon regulatory factor 5

Catalog No: #AP89414



Package Size: #AP89414-1 50ug #AP89414-2 100ug #AP89414-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon regulatory factor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IRF-5
Accession No.	Uniprot ID:Q58DJ0
Uniprot	Q58DJ0
GeneID	615340;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IRF5

Note: This product is for in vitro research use only