

Lymphoid enhancer-binding factor 1

Catalog No: #AP89416



Package Size: #AP89416-1 50ug #AP89416-2 100ug #AP89416-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Lymphoid enhancer-binding factor 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | LEF-1,T cell-specific transcription factor 1-alpha,TCF1-alpha |
| Accession No. | Uniprot ID:Q9UJU2 |
| Uniprot | Q9UJU2 |
| GeneID | 51176; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LEF1

Note: This product is for in vitro research use only