

Endothelial lipase

Catalog No: #AP89418



Package Size: #AP89418-1 50ug #AP89418-2 100ug #AP89418-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelial lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Endothelial cell-derived lipase,EDL,UNQ387/PRO719
Accession No.	Uniprot ID:Q9Y5X9
Uniprot	Q9Y5X9
GeneID	9388;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LIPG

Note: This product is for in vitro research use only