

## Eosinophil lysophospholipase

Catalog No: #AP89421



Package Size: #AP89421-1 50ug #AP89421-2 100ug #AP89421-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eosinophil lysophospholipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Charcot-Leyden crystal protein,CLC,Galectin-10,Gal-10,Lysolecithin acylhydrolase
Accession No.	Uniprot ID:Q05315
Uniprot	Q05315
GeneID	1178;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CLC

Note: This product is for in vitro research use only