

Lymphocyte antigen 6 complex locus protein G6f

Catalog No: #AP89424



Package Size: #AP89424-1 50ug #AP89424-2 100ug #AP89424-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lymphocyte antigen 6 complex locus protein G6f
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G6f
Accession No.	Uniprot ID:Q6MG56
Uniprot	Q6MG56
GeneID	309609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ly6g6f

Note: This product is for in vitro research use only