

Membrane-bound transcription factor site-1 protease

Catalog No: #AP89428



Package Size: #AP89428-1 50ug #AP89428-2 100ug #AP89428-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Membrane-bound transcription factor site-1 protease
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Endopeptidase S1P,Subtilisin/kexin-isozyme 1,SKI-1,KIAA0091,S1P,SKI1
Accession No.	Uniprot ID:Q14703
Uniprot	Q14703
GeneID	8720;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MBTPS1

Note: This product is for in vitro research use only