

## Mitotic spindle assembly checkpoint protein MAD1

Catalog No: #AP89429



Package Size: #AP89429-1 50ug #AP89429-2 100ug #AP89429-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mitotic spindle assembly checkpoint protein MAD1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mitotic arrest deficient 1-like protein 1,MAD1-like protein 1,Mitotic checkpoint MAD1 protein homolog,HsMAD1,Tax-binding protein 181,MAD1,TXBP181
Accession No.	Uniprot ID:Q9Y6D9
Uniprot	Q9Y6D9
GeneID	8379;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MAD1L1

Note: This product is for in vitro research use only