

Stromelysin-2

Catalog No: #AP89432



Package Size: #AP89432-1 50ug #AP89432-2 100ug #AP89432-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Stromelysin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SL-2,Matrix metalloproteinase-10,MMP-10,Transformation-associated protein 34A,Transin-2
Accession No.	Uniprot ID:P07152
Uniprot	P07152
GeneID	117061;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp10

Note: This product is for in vitro research use only