

## Myosin light chain kinase, smooth muscle

Catalog No: #AP89436



Package Size: #AP89436-1 50ug #AP89436-2 100ug #AP89436-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myosin light chain kinase, smooth muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MLCK,Kinase-related protein,KRP,Telokin,MLCK,MLCK1,MYLK1
Accession No.	Uniprot ID:Q15746
Uniprot	Q15746
GeneID	4638;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MYLK

Note: This product is for in vitro research use only