

## Nucleoside diphosphate kinase A

Catalog No: #AP89439



Package Size: #AP89439-1 50ug #AP89439-2 100ug #AP89439-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nucleoside diphosphate kinase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	NDK A,NDP kinase A,nm23-C1,NM23A
Accession No.	Uniprot ID:Q50KA9
Uniprot	Q50KA9
GeneID	480558;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NME1

Note: This product is for in vitro research use only