

OR6B1 Antibody

Catalog No: #44616

Package Size: #44616-1 50ul #44616-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	OR6B1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OR6B1 Antibody detects endogenous levels of total OR6B1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR6B1
Target Name	OR6B1
Other Names	Olfactory receptor 6B1 Olfactory receptor 7-3 OR7-3 Olfactory receptor OR7-9
Accession No.	Swiss-Prot#: O95007NCBI Gene ID:135946
Uniprot	O95007
GeneID	135946;
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only