

## Neprilysin

Catalog No: #AP89441



Package Size: #AP89441-1 50ug #AP89441-2 100ug #AP89441-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neprilysin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Atriopeptidase,Common acute lymphocytic leukemia antigen,CALLA,Enkephalinase,Neutral endopeptidase 24.11,NEP,Skin fibroblast elastase,SFE,CD10,EPN
Accession No.	Uniprot ID:P08473
Uniprot	P08473
GeneID	4311;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MME

Note: This product is for in vitro research use only