

Neuromedin-K receptor

Catalog No: #AP89444



Package Size: #AP89444-1 50ug #AP89444-2 100ug #AP89444-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuromedin-K receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NKR,NK-3 receptor,NK-3R,Neurokinin B receptor,Tachykinin receptor 3,NK3R,TAC3R
Accession No.	Uniprot ID:P29371
Uniprot	P29371
GeneID	6870;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TACR3

Note: This product is for in vitro research use only