

# Bile acid receptor

Catalog No: #AP89445



Package Size: #AP89445-1 50ug #AP89445-2 100ug #AP89445-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bile acid receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Farnesoid X-activated receptor,Farnesol receptor HRR-1,Nuclear receptor subfamily 1 group H member 4
Accession No.	Uniprot ID:Q3SZL0
Uniprot	Q3SZL0
GeneID	540528;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NR1H4

Note: This product is for in vitro research use only