

Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit

Catalog No: #AP89449

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89449-1 50ug #AP89449-2 100ug #AP89449-3 1mg

Description

Product Name	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	DDOST 48 kDa subunit,Oligosaccharyl transferase 48 kDa subunit,OST48
Accession No.	Uniprot ID:Q05052
Uniprot	Q05052
GeneID	404012;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDOST

Note: This product is for in vitro research use only