

Putative phospholipase B-like 2

Catalog No: #AP89457



Package Size: #AP89457-1 50ug #AP89457-2 100ug #AP89457-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Putative phospholipase B-like 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LAMA-like protein 2,Lamina ancestor homolog 2,Phospholipase B domain-containing protein 2
Accession No.	Uniprot ID:Q2KIY5
Uniprot	Q2KIY5
GeneID	514347;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLBD2

Note: This product is for in vitro research use only