

Peroxisome proliferator-activated receptor gamma coactivator 1-alpha

Catalog No: #AP89459

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89459-1 50ug #AP89459-2 100ug #AP89459-3 1mg

Description

Product Name	Peroxisome proliferator-activated receptor gamma coactivator 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PGC-1-alpha,Pgc1,Pgc1a,Ppargc1
Accession No.	Uniprot ID:Q9QYK2
Uniprot	Q9QYK2
GeneID	83516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppargc1a

Note: This product is for in vitro research use only