

Retinoid isomerohydrolase

Catalog No: #AP89470



Package Size: #AP89470-1 50ug #AP89470-2 100ug #AP89470-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Retinoid isomerohydrolase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | All-trans-retinyl-palmitate hydrolase,Retinal pigment epithelium-specific 65 kDa protein,Retinol isomerase |
| Accession No. | Uniprot ID:Q9YGX2 |
| Uniprot | Q9YGX2 |
| GeneID | 395700; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RPE65

Note: This product is for in vitro research use only