

Retinoic acid receptor RXR-beta

Catalog No: #AP89471



Package Size: #AP89471-1 50ug #AP89471-2 100ug #AP89471-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Retinoic acid receptor RXR-beta |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Nuclear receptor coregulator 1,Nuclear receptor subfamily 2 group B member 2,Retinoid X receptor beta,Nr2b2,Rcor-1 |
| Accession No. | Uniprot ID:P49743 |
| Uniprot | P49743 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rxrb

Note: This product is for in vitro research use only