

Protein S100-A9

Catalog No: #AP89472



Package Size: #AP89472-1 50ug #AP89472-2 100ug #AP89472-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	BEE22,Calgranulin-B,Neutrophil cytosolic 23 kDa protein,p23,S100 calcium-binding protein A9
Accession No.	Uniprot ID:P28783
Uniprot	P28783
GeneID	532569;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S100A9

Note: This product is for in vitro research use only