

# Selenoprotein P

Catalog No: #AP89475



Package Size: #AP89475-1 50ug #AP89475-2 100ug #AP89475-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Selenoprotein P
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SeP, Selenoprotein P-like protein
Accession No.	Uniprot ID:P49907
Uniprot	P49907
GeneID	282066;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SELENOP

Note: This product is for in vitro research use only