

## Endophilin-A2

Catalog No: #AP89476

Package Size: #AP89476-1 50ug #AP89476-2 100ug #AP89476-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endophilin-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Endophilin-2,SH3 domain-containing GRB2-like protein 2,SH3p8,SH3P8
Accession No.	Uniprot ID:Q8AXV0
Uniprot	Q8AXV0
GeneID	395202;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SH3GL1

Note: This product is for in vitro research use only